
avantgard

implante dental de conexión cónica
conical connection dental implant

www.ilerimplant.com


Made in Spain

avantgard

implante dental de conexión cónica
conical connection dental implant

Sistema de
implantes
Implant system

Índice - Summary

SISTEMA DE IMPLANTES

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
gmi avantgard conical connection implants are machined in c.p. grade IV titanium of the highest quality, thus ensuring a high chemical stability and excellent biocompatibility. The titanium oxide layer, that is formed on the implant surface after the passivation process, and the careful design of the outer surface of the thread, guarantee good osseointegration and high primary stability.

The external dual inlet thread of the implant helps to reduce surgery time. Moreover, the rounded shapes, self-threading millings and slight tapering facilitate implant insertion, reducing tension at the bone-implant interface and thus preventing threading stress-related problems.

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① Platforms
RP: Standard
WP: Wide

② Lengths

③ Self threading millings
Variable length
3 millings at 120°
Atraumatic tip

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


To promote the adhesion and growth of bone cells, the outer surface of the implant has been treated with the exclusive ADS® treatment (Advanced Double-Grip Surface) which combines a white corundum micro-bubble treatment and acid etching, for a non-smooth roughness that maximizes the contact surface between implant and bone and therefore improves stability and adhesion between them.

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
The ADS® treatment is the result of a series of collaborations with institutes and universities, as well as our own research and tests, that allow us to monitor and validate the optimal biological response to the implant surface: chemical (X-ray scattering studies to find out the composition of the implant and its surface) and topography studies (roughness studies) have been carried out, as well as biological (in vivo response studies in animals) and clinical studies.

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X-ray scattering studies show that the surface of the implants is only made of titanium and oxygen, which are the basic components of the oxide layer that is formed in the implant passivation process and that gives it its excellent corrosion resistance.

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
In order to control and identify the topography of the implant surface, roughness studies were performed with profilometer, as well as with techniques such as SEM (Scanning Electron Microscope) and CLSM (Confocal Laser Scanning Microscope), allowing to obtain and control the mean roughness (Ra) values within the parameters recommended in international publications.




Perfil de rugosidad obtenido con HOMMEL-WERKE T1000 • Roughness profile obtained with HOMMEL-WERKE T1000.

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To determine the biological response of **gmi** implants, in vivo studies were conducted by inserting the implants on animals, leaving them without load during the healing process and performing a histological study to prove the excellent response of the bone cells and the complete osseointegration of the implant.

SEM-BSE: detalles ultraestructurales de las células en biofilm sobre el implante. • SEM-BSE: Ultrastructural details of cells in biofilm on the implant.

Protocolo de inserción de implantes gmi avantgard - Insertion sequence for gmi avantgard implants

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gmi avantgard conical connection implants are packed in a screen-printed cardboard box, with the traceability labels and a double packaging containing the implant system, with the secondary packaging made of technical plastic and the primary packaging made of Pyrex glass to ensure the inert performance of the container. Once packaged in a clean room, they are sterilized with gamma radiation according to regulations.

The new packaging system of the **gmi avantgard** implants consists of a titanium support which houses the implant-implant carrier set, prepared to be inserted directly into the mouth, using a single hex-3.00 mm implant carrier wrench, designed for this purpose. This new system offers the odontologist the following advantages: more control and a better view of the implant insertion process, it is easier to use in reduced interdental spaces and a reduction in surgery time.*

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* Consult the instructions for use in electronic format.

1. Check sticker indicating implant diameter and length and open the tab on the carton.
2. Remove the secondary packaging and product identification labels.
3. Check the integrity of the safety seal and unscrew the plug on the secondary packaging.
4. Remove the primary packaging from the secondary packaging avoiding hitting it against a hard surface.
5. Extract primary plug making a lateral movement, and save it for posterior use.
6. Insert the implant carrier wrench until a slight retention is noticed, addressing the notches on the key with the faces of the implant carrier hex broaching.
7. Check that the key is fully inserted and turn slightly while gently pulling on the assembly.
8. Place the implant in the bone bed by applying a torque of 35 N-cm until the treated area is at the crestal or subcrestal level. Remove the key.
- 9-10. Hold the implant carrier with flat manual wrench and untighten screw using HEX-1.20 mm wrench until the cone is disassembled. Remove the implant carrier.
11. Remove the cover screw from the primary packaging plug, using the HEX-1.20 mm wrench.
12. Manual tightening cover screw to the implant. Remove key and suture the incision.

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

The following is a description of the drilling sequences for the different conical connection implant models, as well as the recommended drilling conditions for their use:

- Lance-shaped: 1200 - 1500 rpm.
- Pilot drills: 900 - 1200 rpm.
- Guide drill: 800 rpm.
- Twist drills:
 - Ø2.80mm → 500 - 700 rpm.
 - Ø3.00 & Ø3.50mm → 400 - 700 rpm.
 - Ø4.00 & Ø5.10mm → 400 - 600 rpm.
- Countersink drills: 200 - 400 rpm.
- For drilling, external cooling with saline solution is required.


avantgard Ø3,75		
	Descr.	Ref.
1	Fresa lanceolada	Lance-shaped drill
2	Fresa piloto	Pilot drill
3	Fresa guía	Guide drill
4	Fresa Ø2,80 mm	Twist drill Ø2.80 mm
5	Fresa Ø3,00 mm	Twist drill Ø3.00 mm
6	Avellanadora Ø3,75 mm	Countersink drill Ø3.75 mm


avantgard Ø4,25		
	Descr.	Ref.
1	Fresa lanceolada	Lance-shaped drill
2	Fresa piloto	Pilot drill
3	Fresa guía	Guide drill
4	Fresa Ø2,80 mm	Twist drill Ø2.80 mm
5	Fresa Ø3,00 mm	Twist drill Ø3.00 mm
6	Fresa Ø3,50 mm	Twist drill Ø3.50 mm
7	Avellanadora Ø4,25 mm	Countersink drill Ø4.25 mm

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



The **gmi avantgard** implant range was designed and tested to be inserted as a single restoration, with a 30-degree maximum inclination to the perpendicular plane to the occlusal plane, in accordance with the following diagram.

		
ø3,75 mm	ø4,25 mm	ø4,75-ø5,75 mm
Incisivos (Excepto centrales maxilares) Incisors (Except central maxillary incisors)	Caninos • Canines Premolares • Premolars Incisivos centrales maxilares Central maxillary incisors	Molares Molars

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The **gmi avantgard** implant range is designed to be inserted applying a 35-40 N-cm torque, although it can bear a maximum tightening torque of 60 N-cm. **gmi** recommends manual insertion of the implant with a dynamometric ratchet wrench to ensure that these values are not exceeded.

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The **gmi avantgard** range of implants and prosthetic attachments are designed and manufactured with tolerances intended to maintain adjustments that minimize the clearance between components and therefore guarantee the stability of the connection. **gmi** recommends the use of original attachments to ensure a perfect fit between the components of the restoration.

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
To prevent prior deformation of the clinic screw, **gmi** recommends restricting its use only for the final tightening of the prosthesis. To make adjustments during the prosthetic procedures, use of the laboratory screw is recommended.

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Once the repair phase of the tissue supporting the implant has finished, in the case of implants fully or partially covered by soft tissue there must be a mucosal channel or tunnel enabling connection of the implant to the prosthesis or secondary structure. Healing abutments are used to form the mucosal route, once the soft tissue has been opened by incision and the cover screw has been withdrawn, by screwing them to the implant.

Key features

- Grade V titanium machining as per ISO-5832 standard.
- Tightening with b/f 1.20 mm hex wrench.
- Available in emergency heights from 1.50 to 4.50 mm.
- Apical bit to facilitate the start of threading.
- Manual tightening torque: Max. 15 N·cm.

Plataforma Platform	Ø cabeza Ø head	Altura (mm) Height (mm)	Ref.
RP	Ø 5,00	1,50	KDPoC4001
		2,50	KDPoC4002
		3,50	KDPoC4003
		4,50	KDPoC4004

Plataforma Platform	Ø cabeza Ø head	Altura (mm) Height (mm)	Ref.
RP	Ø 6,00	2,50	KDPoC4009
		3,50	KDPoC4010
		4,50	KDPoC4011


Plataforma Platform	Ø cabeza Ø head	Altura (mm) Height (mm)	Ref.
WP	Ø 6,00	1,50	KDPoC4005
		2,50	KDPoC4006
		3,50	KDPoC4007
		4,50	KDPoC4008

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When the use of a custom abutment is required, implant carrier abutments can be used. By cutting them in any direction, they can be suited to the needs of the prosthetic. In addition and as a complementary function, implant carrier abutments can be used like impression coping, threading them to the implant screw through the short impression screw and using them to make the negative of the patient's mouth for subsequent positivation on plaster.

Key features

- Grade V titanium machining as per ISO-5832 standard.
- They allow cutting in any direction to suit the needs of the prosthetic.
- Compatible with b/f 3.00 mm hex wrench for use as an implant carrier.
- Available only in non-rotary version.
- Possibility to use as impression coping with short impression screw.
- Premounted clinic screw included.
- Screw tightening torque: 25 N-cm.


Plataforma Platform	Ø hombro Ø shoulder	Ref.
RP	Ø 3,75 - Ø 4,25	KDBoF4001
WP	Ø 4,75	KDBoF4002
	Ø 5,75	KDBoF4003

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When it is necessary to perform a provisional restoration while the bone-implant interfaces or the soft tissues are healing, the temporary abutment is used, defined as attachment that once attached to the crown of the implant with clinic screw, will absorb most of the masticatory stresses while avoiding overloading the implant while the osseointegration process is completed.

Key features

- Machined in high strength technical plastic.
- Rotary version only for multiple restorations.
- Anti-rotational plane and transverse retention to provide good fixation.
- Clinic screw included.
- Manual tightening torque: Max. 15 N-cm.

Plataforma Platform	Tipo conexión Connection type	Ref.
RP	Anti-rotatoria • Non-rotary	KDHoF4001
	Rotatoria • Rotary	KDHoF4002
WP	Anti-rotatoria • Non-rotary	KDHoF4003
	Rotatoria • Rotary	KDHoF4004

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When it is necessary to correct an extreme angular position of the implant inside the bone structure in relation to adjacent natural or prosthetic teeth, angled abutments should be used. Once mounted and oriented on the implant, these abutments will allow the parallelization of the new crown with other structures.

Key features

- Grade V titanium machining as per ISO-5832 standard.
- Allow non-parallelism correction up to 20 degrees.
- Clinic screw included.
- Abutment shoulder adapted to the shape of the gum.
- Screw tightening torque: 25 N·cm.

Plataforma Platform	Ang. (°)	Altura (mm) Height (mm)	Ref.
RP	15°	2,50	KDMoF4013
	15°	4,50	KDMoF4014
	20°	2,50	KDMoF4015
	20°	4,50	KDMoF4016

Plataforma Platform	Ang. (°)	Altura (mm) Height (mm)	Ref.
WP	15°	2,50	KDMoF4017
	15°	4,50	KDMoF4018
	20°	2,50	KDMoF4019
	20°	4,50	KDMoF4020

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

Once the mucosal tunnel has been formed by the healing abutment and after removing it from this location, it is appropriate to use the abutment, which is defined as the machined tapered attachment that, after being fixed with the clinic screw to the top of the implant, acts as a support for the clinic crown and is attached to it using cement.

Key features

- Grade V titanium machining as per ISO-5832 standard.
- Anti-rotational planes to increase retention.
- Available only in non-rotary version.
- Clinic screw included.
- Available in four emergency heights from 1.50 to 4.50 mm.
- Screw tightening torque: 25 N·cm.

Plataforma Platform	Ø hombro Ø shoulder	Altura (mm) Height (mm)	Ref.
RP	Ø 4,00	1,50	KDMoF4001
		2,50	KDMoF4002
		3,50	KDMoF4003
		4,50	KDMoF4004
	Ø 5,00	1,50	KDMoF4005
		2,50	KDMoF4006
		3,50	KDMoF4007
		4,50	KDMoF4008

Plataforma Platform	Ø hombro Ø shoulder	Altura (mm) Height (mm)	Ref.
WP	Ø 6,00	1,50	KDMoF4009
		2,50	KDMoF4010
		3,50	KDMoF4011
		4,50	KDMoF4012

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
Once the mucosal tunnel has been formed by the healing abutment and after removing the cap from this location, it is appropriate to use the titanium base, which is defined as the machined tapered attachment that, after being fixed with the clinic screw to the top of the implant, acts as a support for a metallic or ceramic meso-structure casted or machined with CAD/CAM.

Key features



- Grade V titanium machining as per ISO-5832 standard.
- Non-rotary and rotary version.
- Clinic screw included.
- Available in emergency height of 1.50 mm and 3.00 mm.
- Non-rotary titanium base with indexation.
- Screw tightening torque: 25 N-cm.
- Scan directly on Titanium base or alternative scan on Scan abutment.

Base titanio - Titanium base			
Plataforma Platform	Altura (mm) Height (mm)	Tipo conexión Connection type	Ref.
RP	1,50	Anti-rotatoria • Non-rotary	KDMoF4021
	1,50	Rotatoria • Rotary	KDMoF4023
	3,00	Anti-rotatoria • Non-rotary	KDMoF4025
	3,00	Rotatoria • Rotary	KDMoF4027
WP	1,50	Anti-rotatoria • Non-rotary	KDMoF4022
	1,50	Rotatoria • Rotary	KDMoF4024
	3,00	Anti-rotatoria • Non-rotary	KDMoF4026
	3,00	Rotatoria • Rotary	KDMoF4028

Calcinable - Castable		
Plataforma Platform	Altura total (mm) Total Height (mm)	Ref.
RP	11,00	KDCoF3605
WP	11,00	KDCoF3606

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Following the first stage of surgery and therefore after implant osseointegration is achieved, the multi-esthetic abutment (MEA) is used, that is defined as the machined attachment that, after being attached to the crown of the implant, serves as a support for a multiple restoration, thus decreasing the height between the implant connection and the structure connection, and providing a passive adjustment of the structure.

Key features

- Grade V titanium machining as per ISO-5832 standard.
- Available in emergency heights from 1.50 to 4.50 mm.
- One-piece straight abutments in rotary version.
- Tightened with hex wrench: B/F 2.00 mm (KYLoCo149).
- Impression coping with screw included.
- Castable with clinic screw included.
- Temporary cylinder in grade V titanium with clinic screw included.
- Angled multi-esthetic abutments allow non-parallelism up to 30°.
- Tightening torques: ME straight abut. (30 N-cm), ME angled abut. (20 N-cm) and healing caps (Manual: Max. 15 N-cm).

Pilares multi-estético gmi avantgard - gmi avantgard multi-esthetic abutments

Pilares multi-estético inclinados Multi-esthetic angled abutments

Plataforma Platform	Ang. (°)	Altura (mm) Height (mm)	Ref.
RP	17°	2,50	KDGoF4009
	17°	3,50	KDGoF4010
	30°	4,00	KDGoF4011

Pilares multi-estético rectos Multi-esthetic straight abutments

Plataforma Platform	Altura (mm) Height (mm)	Ref.
RP	1,50	KDGoF4001
	2,50	KDGoF4002
	3,50	KDGoF4003
	4,50	KDGoF4004
WP	1,50	KDGoF4005
	2,50	KDGoF4006
	3,50	KDGoF4007
	4,50	KDGoF4008

Pilares cicatrización PME Healing abutments MEA

Plataforma Platform	Ø Diámetro (mm) Ø Diameter (mm)	Ref.
RP - WP	Ø 4,80	KDPoF0035
	Ø 6,00	KDPoF0037


Tomas impresión cubeta abierta (CA) larga PME Open tray (OT) long impression coping MEA

Plataforma Platform	Ref.
RP - WP	KDSOf0015

Pilares multi-estético gmi avantgard - gmi avantgard multi-esthetic abutments

Tomas impresión cubeta cerrada (CC) PME
Closed tray (CT) impression coping MEA

Plataforma Platform	Cantidad (Ud.) Quantity (Un.)	Ref.
RP - WP	1	KDS0F0017
	5	KDS0F5001

Cilindro provisional PME
Temporary cylinder MEA

Plataforma Platform	Ref.
RP - WP	KDGoF0028


Calcinable PME
Castable MEA

Plataforma Platform	Ref.
RP - WP	KDCoF0117

Base titanio PME
Titanium base MEA

Plataforma Platform	Ref.
RP - WP	KDGoF0035

Scanbody de titanio PME
Titanium scanbody MEA

Plataforma Platform	Ref.
RP - WP	KDloF0002

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
Once the insertion of the implant in the bone bed has been completed, the bioesthetic base is used, which is defined as the machined attachment that, once fixed to the coronal part of the implant, allows single restorations at the abutment level.

Placing the bioesthetic base during surgery allows the soft tissues to heal around it.

Once the tissue has healed, it serves as the basis for the rehabilitation, avoiding the connection and disconnection of the components to the implant and, therefore, keeping the biological seal intact.

Key features

- Grade V titanium machining as per ISO-5832 standard.
- Available in emergency heights from 2.50 to 4.50 mm.
- Straight two-piece bases, including the screw.
- Tightened with hex wrench: B/F 2.00 mm (KYLoC0149).
- Impression coping with screw included.
- Castable.
- Tightening torques: BE base (25 N-cm) and healing abutment (Manual: Max. 15 N-cm).

Bases bioestéticas Bioesthetic bases		
Plataforma Platform	Altura (mm) Height (mm)	Ref.
RP	2,50	KDGoF4012
	3,50	KDGoF4013
	4,50	KDGoF4014
WP	2,50	KDGoF4015
	3,50	KDGoF4016
	4,50	KDGoF4017

Pilar cicatrización BE BE healing abutment		
Plataforma Platform	Ø Diámetro (mm) Ø Diameter (mm)	Ref.
RP - WP	Ø 4,80	KDPoC0038

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



Once the transfer of the position has been made using the impression coping from the biological medium to the laboratory working model, it is necessary to use replicas, which are defined as the machined attachments designed on one hand to reproduce the implant position in the mouth of the patient to a working model, and on the other to serve as a connection model for the construction in the laboratory of the structure to carry out the restoration.

Key features

- Grade V titanium machining as per ISO-5832 standard.
- Anti-rotation millings and vertical retention neck for a good fixation.

Réplicas de implante Implant replicas

Plataforma Platform	Ref.
RP	KDRoC4001
WP	KDRoC4002

Réplicas digitales Digital replicas

Plataforma Platform	Ref.
RP	KDRoF4003
WP	KDRoF4004

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Plataforma Platform	Ref.
RP - WP	KDRoC0012

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Plataforma Platform	Ref.
RP - WP	KDRoF0014

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
Impression coping are defined as the machined attachments that, connecting to the coronal part of the implant within the oral cavity and axially fixed by the impression screw, is used to make the transfer of the position of the implant in the biological medium to a laboratory model, on which the prosthetic will work to create the dental crowns or prostheses. This process is achieved using impression materials placed in an appropriate manner, that will harden in the oral cavity and produce the negative of the patient's mouth. Once extracted, replicas are placed and the casting plaster is poured to obtain the positive mould where the replica will be placed in the original position of the implant in the mouth.

Key features

- Machined open tray impression coping in grade V titanium according to ISO-5832.
- Impression screws in grade V titanium according to ISO-5832.
- Tightening of impression screws with b/f 1.20 mm hex wrench.
- Apical bit to facilitate the start of threading impression screw.
- Anti-rotation millings and vertical retention neck for a good fixation.
- Impression coping screw included.
- Long impression screw with longitudinal millings for manual tightening.
- CT impression coping machined in technical plastic.
- Manual tightening torque: Max. 15 N-cm.

Tomas impresión cubeta abierta (CA) Open tray (OT) impression coping		
Plataforma Platform	Tornillo Screw	Ref.
RP	Largo • Long	KDSOf4001
	Corto • Short	KDSOf4003
WP	Largo • Long	KDSOf4002
	Corto • Short	KDSOf4004

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Tomas impresión cubeta cerrada (CC) Closed tray (CT) impression coping

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Key features

- Machined grade V titanium screws according to ISO-5832.
- Tightening of screws with b/f 1.20 mm hex wrench.

Tornillo clínica Clinic screw

Plataforma Platform	Hex. llave (mm) Hex. wrench (mm)	Ref.
RP	⊕ 1,20 mm	KDToC4003
WP	⊕ 1,20 mm	KDToC4004

Tornillo laboratorio Laboratory screw

Plataforma Platform	Hex. llave (mm) Hex. wrench (mm)	Ref.
RP	⊕ 1,20 mm	KDToC4005
WP	⊕ 1,20 mm	KDToC4006

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MEA and BE clinical screws are defined as machined devices used for the final anchoring of the ME and BE abutments, while laboratory screws are the machined devices used for the temporary anchoring of ME and BE abutments to the replica during the fabrication of the prostheses in the laboratory.

Key features

- Machined grade V titanium screws according to ISO-5832.
- Tightening of screws with b/f 1.20 mm and 2.00 mm hex wrench.

Tornillo PME MEA screw	
Plataforma Platform	Ref.

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
gmi avantgard surgical kit consists of an autoclavable case made of high strength resin, containing all the necessary material for the correct positioning of **gmi avantgard** conical connection dental implant system. The components are arranged in such a way to be easy to find in each moment and in any stage of the intervention.

gmi surgical kits are distributed in cartons, made of high strength technical polymer, and should be steam autoclaved before use at a maximum temperature of 134 °C.

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
In the chart are detailed all the symbols appearing on the **gmi avantgard** implant system's labelling and packaging and their corresponding description.


Descripción de los símbolos del sistema de implantes avantgard Description of symbols of avantgard implants system

- | | | |
|----|--|---|
| 1 | REF | Referencia - Reference |
| 2 | LOT | Número de Lote - Batch code |
| 3 |  | Fecha de caducidad - Used-by-date |
| 4 |  | Fabricante - Manufacturer |
| 5 | 
2797 | Marco de conformidad europeo
European conformity mark |
| 6 | STERILE | Esterilizado usando radiación gamma
Sterilized using gamma radiation |
| 7 | 
ilerimplant.com/IFU | Instrucciones de uso electrónicas
Electronic instructions for use |
| 8 |  | No utilizar si el envase está dañado o abierto
Do not use if packaging is damaged or opened |
| 9 |  | De un solo uso - Single use only |
| 10 | 
NON STERILE | No estéril - Non-sterile |
| 11 | Rx Only | Venta limitada a dentistas o médicos autorizados
Sale by or on the order of a licensed dentist |

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