
frontier

implante dental de conexión interna hexagonal
dental implant of hexagonal internal connection

www.ilerimplant.com


Made in Spain

frontier

implante dental de conexión interna hexagonal
dental implant of hexagonal internal connection

Sistema de
implantes
Implant system

Índice - Summary

SISTEMA DE IMPLANTES

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
gmi frontier internal connection hexagonal implants are machined in c.p. grade IV titanium of the highest quality, thus ensuring a high chemical stability and excellent biocompatibility. The titanium oxide layer, that is formed on the implant surface after the passivation process, and the careful design of the outer surface of the thread, guarantee good osseointegration and high stability.

The external dual inlet thread of the implant helps to reduce surgery time. Moreover, the rounded shapes, self-threading millings and slight tapering facilitate implant insertion, reducing tension at the bone-implant interface and thus preventing threading stress-related problems.

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

To promote the adhesion and growth of bone cells, the outer surface of the implant has been treated with the exclusive ADS® treatment (Advanced Double-Grip Surface) which combines a white corundum micro-bubble treatment and acid etching, for a non-smooth roughness that maximizes the contact surface between implant and bone and therefore improves stability and adhesion between them.

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
The ADS® treatment is the result of a series of collaborations with institutes and universities, as well as our own research and tests, that allow us to monitor and validate the optimal biological response to the implant surface: chemical (X-ray scattering studies to find out the composition of the implant and its surface) and topography studies (roughness studies) have been carried out, as well as biological (in vivo response studies in animals) and clinical studies.

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X-ray scattering studies show that the surface of the implants is only made of titanium and oxygen, which are the basic components of the oxide layer that is formed in the implant passivation process and that gives it its excellent corrosion resistance.

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
In order to control and identify the topography of the implant surface, roughness studies were performed with profilometer, as well as with techniques such as SEM (Scanning Electron Microscope) and CLSM (Confocal Laser Scanning Microscope), allowing to obtain and control the mean roughness (Ra) values within the parameters recommended in international publications.




Perfil de rugosidad obtenido con HOMMEL-WERKE T1000 • Roughness profile obtained with HOMMEL-WERKE T1000.

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To determine the biological response of **gmi** implants, in vivo studies were conducted by inserting the implants on animals, leaving them without load during the healing process and performing a histological study to prove the excellent response of the bone cells and the complete osseointegration of the implant.

SEM-BSE: detalles ultraestructurales de las células en biofilm sobre el implante. • SEM-BSE: Ultrastructural details of cells in biofilm on the implant.

Protocolo de inserción de implantes gmi frontier - Insertion sequence for gmi frontier implants

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gmi frontier internal connection implants are packed in a screen-printed cardboard box, with the traceability labels and a double packaging containing the implant system, with the secondary packaging made of technical plastic and the primary packaging made of Pyrex glass to ensure the inert performance of the container. Once packaged in a clean room, they are sterilized with gamma radiation according to regulations.

The new packaging system of the **gmi frontier** implants consists of a titanium support which only houses the dental implant, prepared to be inserted directly into the mouth, using a single hex-2.45 mm implant carrier wrench, designed for this purpose. This new system offers the odontologist the following advantages: more control and a better view of the implant insertion process, it is easier to use in reduced interdental spaces and a reduction in surgery time.*

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
* Consult the instructions for use in electronic format.


1. Check sticker indicating implant diameter and length and open the tab on the carton.
2. Remove the secondary packaging and product identification labels.
3. Check the integrity of the safety seal and unscrew the plug on the secondary packaging.
4. Remove the primary packaging from the secondary packaging avoiding hitting it against a hard surface.
5. Extract and save the primary packaging plug making a lateral movement.
6. Insert the HEX-2.45 mm implant carrier wrench until a slight retention is noticed, addressing the hex of the key with the faces of the implant hex broaching.
7. Check that the key is fully inserted and turn slightly while gently pulling on the assembly.
8. Place the implant in the bone bed by applying a torque of 35 N-cm until the treated area is at the crestal or subcrestal level. Remove the key.
9. Remove the cover screw from the primary packaging plug, using the HEX-1.20 mm wrench.
10. Screw the cover screw to the implant by applying a maximum torque of 15 N-cm. Remove key and suture the incision.

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
The following is a description of the drilling sequences for the different internal connection implant models, as well as the recommended drilling conditions for their use:


- Lance-shaped: 1200 - 1500 rpm.
- Pilot drills: 900 - 1200 rpm.
- Guide drill: 800 rpm.
- Twist drills:
 - Ø2.80mm → 500 - 700 rpm.
 - Ø3.00 & Ø3.50mm → 400 - 700 rpm.
 - Ø4.00 & Ø5.10mm → 400 - 600 rpm.
- Countersink drills: 200 - 400 rpm.
- For drilling, external cooling with saline solution is required.

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






The **gmi frontier** implant range was designed and tested to be inserted as a single restoration, with a 30-degree maximum inclination to the perpendicular plane to the occlusal plane, in accordance with the following diagram.

		
ø3,30-ø3,75 mm	ø4,25 mm	ø4,75-ø5,75 mm
Incisivos (Excepto centrales maxilares) Incisors (Except central maxillary incisors)	Caninos • Canines Premolares • Premolars Incisivos centrales maxilares Central maxillary incisors	Molares Molars

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The **gmi frontier** implant range is designed to be inserted applying a 35-40 N-cm torque, although it can bear a maximum tightening torque of 60 N-cm. **gmi** recommends manual insertion of the implant with a dynamometric ratchet wrench to ensure that these values are not exceeded.

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The **gmi frontier** range of implants and prosthetic attachments are designed and manufactured with tolerances intended to maintain adjustments that minimize the clearance between components and therefore guarantee the stability of the connection. **gmi** recommends the use of original attachments to ensure a perfect fit between the components of the restoration.

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
To prevent prior deformation of the clinic screw, **gmi** recommends restricting its use only for the final tightening of the prosthesis. To make adjustments during the prosthetic procedures, use of the laboratory screw is recommended.

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Once the repair phase of the tissue supporting the implant has finished, in the case of implants fully or partially covered by soft tissue there must be a mucosal channel or tunnel enabling connection of the implant to the prosthesis or secondary structure. Healing abutments are used to form the mucosal route, once the soft tissue has been opened by incision and the cover screw has been withdrawn, by screwing them to the implant.

Key features

- Grade V titanium machining as per ISO-5832 standard.
- Tightening with b/f 1.20 mm hex wrench.
- Alternative tightening with a flat tip screwdriver.
- Available in emergency heights from 2.00 to 6.00 mm.
- Apical bit to facilitate the start of threading.

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
When the use of a custom abutment is required, millable abutments can be used. By cutting them in any direction, they can be suited to the needs of the prosthetic. In addition and as a complementary function, millable implant abutments can be used like impression coping, threading them to the implant screw through the short impression screw and using them to make the negative of the patient's mouth for subsequent positivation on plaster.

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Key features

- Grade V titanium machining as per ISO-5832 standard.
- They allow cutting in any direction to suit the needs of the prosthetic.
- Compatible with b/f 4.00 mm hex wrench for use as an implant carrier.
- Available only in non-rotary version.
- Possibility to use as impression coping with short impression screw.
- Clinic screw part no.: KDToc3601 included.


Plataforma Platform	Ref.
RP	KDBoF3601
WP	KDBoF3602

Pilares inclinados estéticos gmi frontier - gmi frontier esthetic angled abutments

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When it is necessary to correct an extreme angular position of the implant inside the bone structure in relation to adjacent natural or prosthetic teeth, angled abutments should be used. Once mounted and oriented on the implant, these abutments will allow the parallelization of the new crown with other structures.

Key features

- Grade V titanium machining as per ISO-5832 standard.
- Allow non-parallelism correction up to 20 degrees.
- Clinic screw part no.: KDT0C3601 included.
- Abutment shoulder adapted to the shape of the gum.

Pilares inclinados estéticos Esthetic angled abutments		
Plataforma Platform	Ang. (°)	Ref.
RP	15°	KDMoF3607
	20°	KDMoF3611
WP	15°	KDMoF3609
	20°	KDMoF3612

Pilares inclinados Angled abutments		
Plataforma Platform	Ang. (°)	Ref.
RP	20°	KDMoF3613
WP	20°	KDMoF3614

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Once the mucosal tunnel has been formed by the healing abutment and after removing it from this location, it is appropriate to use the abutment, which is defined as the machined tapered attachment that, after being fixed with the clinic screw to the top of the implant, acts as a support for the clinic crown and is attached to it using cement.


Key features

- Grade V titanium machining as per ISO-5832 standard.
- Anti-rotational planes to increase retention.
- Abutments in b/f 2.45 mm hex non-rotary version.
- Clinic screw part no.: KDT0C3601 included.
- Available in three emergency heights from 1.50 to 3.50 mm.

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
When it is necessary to perform a provisional restoration while the bone-implant interfaces or the soft tissues are healing, the temporary cylinder is used, defined as technical plastic machined attachment that once attached to the crown of the implant with clinic screw, will absorb most of the masticatory stresses while avoiding overloading the implant while the osseointegration process is completed.

Key features

- Machined in high strength technical plastic.
- Abutment shoulder adapted to the shape of the gum.
- Anti-rotational plane and transverse retention to provide good fixation.
- Clinic screw part no.: KDToc3601 included.

Plataforma Platform	Tipo conexión Connection type	Ref.
RP	Anti-rotatoria • Non-rotary	KDHOf3601
	Rotatoria • Rotary	KDHOf3603
WP	Anti-rotatoria • Non-rotary	KDHOf3602
	Rotatoria • Rotary	KDHOf3604

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
Once the mucosal tunnel has been formed by the healing abutment and after removing the cap from this location, it is appropriate to use the titanium base, which is defined as the machined tapered attachment that, after being fixed with the clinic screw to the top of the implant, acts as a support for a metallic or ceramic meso-structure casted or machined with CAD/CAM.

Key features

- Grade V titanium machining as per ISO-5832 standard.
- Non-rotary and rotary version.
- Clinic screw part no.: KDT0C3601 included.
- Available in emergency height of 0.60 mm
- Non-rotary titanium base with indexation.
- Scan directly on Titanium base or alternative scan on Scan abutment.

Base titanio - Titanium base

Plataforma Platform	Altura (mm) Height (mm)	Tipo conexión Connection type	Ref.
RP	0,65	Anti-rotatoria • Non-rotary	KDMoF3615
	0,65	Rotatoria • Rotary	KDMoF3616
	2,00	Anti-rotatoria • Non-rotary	KDMoF3619
	2,00	Rotatoria • Rotary	KDMoF3620
WP	0,65	Anti-rotatoria • Non-rotary	KDMoF3617
	0,65	Rotatoria • Rotary	KDMoF3618
	2,00	Anti-rotatoria • Non-rotary	KDMoF3621
	2,00	Rotatoria • Rotary	KDMoF3622


Calcinable - Castable

Plataforma Platform	Altura total (mm) Total Height (mm)	Ref.
RP	11,00	KDCoF3605
WP	11,00	KDCoF3606

Pilar escaneado - Scanabutment

Plataforma Platform	Altura (mm) Height (mm)	Ref.
RP	5,30	KDloF3601
WP	5,30	KDloF3602

Base titanio gmi frontier - gmi frontier titanium base

Cuerpo pilar escaneado - Scanbody

Plataforma Platform	Altura (mm) Height (mm)	Ref.
RP - WP	7,50	KDIOc4003

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Once the healing abutment has formed the mucosal tunnel and the healing abutment has been removed from this location, the corresponding ball abutment is used, which is defined as a machined attachment being fixed to the top of the implant and serving as the retention element for the removable prosthesis, with the previous addition of the retention caps that will be fixed to the top spherical part of the ball abutment.

Key features

- Grade V titanium machining as per ISO-5832 standard.
- Available in emergency heights from 1.00 to 4.00 mm.
- Pink retention cap and metal container included.
- Protection disk included.
- Tightening with b/f 2.45 mm hex wrench, part no.: KYLoC0065.
- Instruments for caps: Insertion / removal tool for caps (RH485IC).

Plataforma Platform	Altura (mm) Height (mm)	Ref.
RP	1,00	KDFoF3609
	2,00	KDFoF3610
	3,00	KDFoF3611
	4,00	KDFoF3612

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Following the first stage of surgery and therefore after implant osseointegration is achieved, the multi-esthetic abutment (MEA) is used, that is defined as the machined attachment that, after being attached to the crown of the implant, serves as a support for a multiple restoration, thus decreasing the height between the implant connection and the structure connection, and providing a passive adjustment of the structure.

Key features

- Grade V titanium machining as per ISO-5832 standard.
- Available in emergency heights from 1.00 to 4.00 mm.
- One-piece straight abutments in rotary model.
- Tightened with hex keys:
 - RP: B/F 2.00 mm (KYLoC0149)
 - WP: B/F 2.50 mm (KYLoC0098)
- Impression coping with screw included.
- Castable with clinic screw included.
- Temporary cylinder in grade V titanium with clinic screw included.
- Angled multi-esthetic abutments allow non-parallelism up to 30°.

Pilares multi-estético gmi frontier - gmi frontier multi-esthetic abutments

Pilares multi-estético inclinados Multi-esthetic angled abutments

Plataforma Platform	Ang. (°)	Altura (mm) Height (mm)	Ref.
RP	17°	2,50	KD6oF3609
	17°	3,50	KD6oF3610
	30°	4,00	KD6oF3611

Pilares multi-estético rectos Multi-esthetic straight abutments

Plataforma Platform	Altura (mm) Height (mm)	Ref.
RP	1,00	KD6oC3601
	2,00	KD6oC3602
	3,00	KD6oC3603
	4,00	KD6oC3604
WP	1,00	KD6oC3605
	2,00	KD6oC3606
	3,00	KD6oC3607
	4,00	KD6oC3608


Pilares cicatrización PME Healing abutments MEA

Plataforma Platform	Ø Diámetro (mm) Ø Diameter (mm)	Ref.
RP	Ø 4,80	KDPoF0035
	Ø 6,00	KDPoF0037
WP	Ø 6,00	KFPoF0036

Tomas impresión cubeta abierta (CA) larga PME Open tray (OT) long impression coping MEA

Plataforma Platform	Ref.
RP	KDSoF0015
WP	KDSoF0016

Pilares multi-estético gmi frontier - gmi frontier multi-esthetic abutments

Toma impresión cubeta cerrada (CC) PME Closed tray (CT) impression coping MEA

Plataforma Platform	Cantidad (Ud.) Quantity (Un.)	Ref.
RP	1	KDSOf0017
	5	KDSOf5001

Cilindros provisionales PME Temporary cylinder MEA

Plataforma Platform	Ref.
RP	KDGoF0028
WP	KDGoF0029


Calcinables PME Castables MEA

Plataforma Platform	Ref.
RP	KDCoF0117
WP	KDCoF0118

Base titanio PME Titanium base MEA

Plataforma Platform	Ref.
RP	KDGoF0035
WP	KDGoF0036

Scanbody de titanio PME Titanium scanbody MEA

Plataforma Platform	Ref.
RP	KDioF0002
WP	KDioF0003

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



Once the transfer of the position has been made using the impression coping from the biological medium to the laboratory working model, it is necessary to use replicas, which are defined as the machined attachments designed on one hand to reproduce the implant position in the mouth of the patient to a working model, and on the other to serve as a connection model for the construction in the laboratory of the structure to carry out the restoration.

Key features

- Machined replicas in high quality stainless steel.
- Anti-rotation millings and vertical retention neck for a good fixation.

Réplicas de implante Implant replicas	
Plataforma Platform	Ref.
RP	KDRoC3601

WP	KDRoC3602
----	-----------

Réplicas digitales Digital replicas	
Plataforma Platform	Ref.
RP	KDRoF3603

WP	KDRoF3604
----	-----------

Réplicas para pilares multi-estético MEA replicas	
Plataforma Platform	Ref.
RP	KDRoC0012

WP	KDRoC0013
----	-----------

Réplicas digitales para pilares multi-estético MEA digital replicas	
Plataforma Platform	Ref.
RP	KDRoF0014

WP	KDRoF0015
----	-----------

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Impression coping are defined as the machined attachments that, connecting to the coronal part of the implant within the oral cavity and axially fixed by the impression screw, are used to make the transfer of the position of the implant in the biological medium to a laboratory model, on which the prosthetic will work to create the dental crowns or prostheses. This process is achieved using impression materials placed in an appropriate manner, that will harden in the oral cavity and produce the negative of the patient's mouth. Once extracted, replicas are placed and the casting plaster is poured to obtain the positive mould where the replica will be placed in the original position of the implant in the mouth.

Key features


- Machined open tray impression coping in grade V titanium according to ISO-5832.
- Impression screws in grade V titanium according to ISO-5832.
- Tightening of impression screws with b/f 1.20 mm hex wrench.
- Apical bit to facilitate the start of threading impression screw.
- Anti-rotation millings and vertical retention neck for a good fixation.
- Impression coping screw included.
- Long impression screw with longitudinal millings for manual tightening.
- CT impression coping machined in technical plastic.
- CT impression coping whole formed by CT impression coping + hexagonal abutment + screw.

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Tomas impresión cubeta cerrada (CC) Closed tray (CT) impression coping

Plataforma Platform	Ref.
RP	KDSOf3606
WP	KDSOf3607
Cantidad (Ud.) Quantity (Un.)	Ref.
5	KDSOf0014

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
To replace a dental crown with a restoration dental prosthesis, a castable is used, which is defined as a meltable plastic machined attachment that, once positioned over the coronal part of the internal connection implant, will be the initial basis to model the form of the final structure, to be cast in the corresponding metal.

Key features

- Machined castables in white meltable plastic material.
- Physiologically inert and excellent dimensional stability.
- Available in rotary and non-rotary version.
- Clinic screw part no.: KDTOC3601 included.

Plataforma Platform	Tipo de conexión Connection type	Ref.
RP	Anti-rotatoria • Non-rotary	KDCOf3611
	Rotatoria • Rotary	KDCOf3612
WP	Anti-rotatoria • Non-rotary	KDCOf3613
	Rotatoria • Rotary	KDCOf3614

Tornillos clínica y laboratorio gmi frontier - gmi frontier clinic and laboratory screws

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Clinic screws are defined as machined devices used for the final anchoring of the abutments, the millable implant abutment or the castable directly to the implant, while laboratory screws are the machined devices used for the temporary anchoring of castables to the replica of the implant used by the dental laboratory where the construction of dental crowns or prostheses takes place.

Key features

- Machined grade V titanium screws according to ISO-5832.
- Tightening of screws with b/f 1.20 mm hex wrench.
- Apical bit to facilitate the start of threading.

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gmi frontier surgical kit consists of an autoclavable case made of high strength resin, containing all the necessary material for the correct positioning of gmi frontier internal connection dental implant system. The components are arranged in such a way to be easy to find in each moment and in any stage of the intervention.

gmi surgical kits are distributed in cartons, made of high strength technical polymer, and should be steam autoclaved before use at a maximum temperature of 134°C.

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In the chart are detailed all the symbols appearing on the **gmi frontier** implant system's labelling and packaging and their corresponding description.

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